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Instructions on the individual manufacture of orthopaedic inlays using the HAIX® CO System

(HAIX® Certified Orthopedic System)

Only a few points need to be borne in mind with regard to the manufacture of HAIX® inlay soles in compliance with Rule 112-191 of the German Association of Statutory Accident Insurance Providers (DGUV) and in accordance with Austrian Standard ÖN-Z1259-2012.

The HAIX® CO System comprises a robust green foam material with a resistant cover bearing the HAIX® Logo Design. This prefabricated foam can be used for shoe sizes up to EU 51 / UK 15.

You can build up and adjust the inlay to meet the requirements of the wearer as long as you take account of the area beneath the toe cap which cannot be altered (at least 10 centimetres from the tip of the footbed to the area of the balls of the feet).

PLEASE NOTE:

Only an EVA wedge with a minimum hardness of 40° shore A and a maximum hardness of 50° shore A may be used to build up the footbed.

OUR RECOMMENDATION:

1. EVA blank approximately 40° shore A –
for safety and work shoes subjected to normal use
2. EVA blank approximately 50° shore A –
for safety and work shoes subjected to heavy-duty use (e.g. fire fighter with equipment)

WARNING:

The antistatic properties and ESD capability of all EG prototype-tested HAIX® professional and safety footwear are secured via an electrically conductible textile strip on the inside of the shoe. This grey anti-static band must extend at least 5 centimetres into the shoe. In order to continue to be able to guarantee conductivity, the inlay must not cover the strip and the contact area to the foot needs to be about 4 centimetres. In light of the above information regarding the conductive strip, the use of antistatic or ESD-capable materials (such as EVA/adhesive) is NOT NECESSARY.

The authorization applies to the HAIX shoes, which are marked with the CO-System.

"HAIX® shoes without CE marking are exempted from compliance with the German Social Accident Insurance (DGUV) rule 112-191 (BGR 191) and application of the ÖN-Z1259-2012."

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Hermann Springer GmbH
COMFORT A/S